

# Shred-n'Rich™

PORTION-CONTROLLED NESTING & ENRICHMENT MATERIAL

## NESTING + ENRICHMENT AIDS IN THERMOREGULATION AND DECREASES ANIMAL STRESS

### PORTION CONTROLLED AND EASY TO DISPENSE

Ensures portion-controlled distribution of enrichment material, helping to keep per diem costs on track. Material may be transferable at cage change, providing habitat acclimation for animals and further reducing overall costs. With no pulling, tearing or prepping required, Shred-n'Rich can be removed from packaging and immediately dispensed into cages.

### SUPERIOR ABSORBENCY

Absorbs seven times its own weight of aqueous liquids

### ENHANCED CAGE ENVIRONMENT

Shred-n'Rich provides the perfect environmental, physical and mental stimulant for rodents and burrowing animals.

### BRIGHT WHITE AND DUST-FREE

The bright white color aids in post-operative or toxicology observations. As a virtually dust-free product, Shred-n'Rich is perfect for use with nude mice.

### AVAILABLE OPTIONS

Shred-n'Rich is available irradiated or non-irradiated. Certified product and its documentation is available as a value added option. This product can be autoclaved in-situ or in bulk.

### PRODUCT PURITY

Shred-n'Rich is made from a highly absorbent material that is used in a variety of FDA regulated industries including Food Packaging.

### PRODUCT SPECIFICATIONS AND PACKAGING

Shred-n'Rich is sold by the case.

Product Dimensions	2" x 4" x 0.143"
Eaches per case	3000
Case size	23.25" x 11" x 19"
Case weight (packaging & contents)	11.23 lbs.

# PROVIDING A COMPLETE LINE OF BEDDING, NESTING AND ENRICHMENT PRODUCTS

	PRODUCT	DESCRIPTION	PRODUCT CODE	CAGE FILLS	PACKAGING	PALLET STACK
BEDDING	<b>Bed-o'Cobs® 1/4"</b>	1/4" Cob	4B	115	1.25 cu ft/35 L bag	70
	<b>Bed-o' Cobs® 1/8"</b>	1/8" Cob	8B	115	1.25 cu ft/35 L bag	60
	<b>Bed-o' Cobs® Combo</b>	1/8" & 1/4 Cob	CB	115	1.25 cu ft/35 L bag	70
	<b>Pure-o' Cel®</b>	Paper Squares	PB	115	1.25 cu ft/35 L bag	80
	<b>XTRAFine™</b>	Finely Ground Extra Absorbent Cob	XFB	115	1.25 cu ft/35 L bag	60
	<b>Pel-o'Cobs®</b>	Pelleted Cob Bedding	PELB	-	1.25 cu ft/35 L bag	50
BEDDING + ENRICHMENT	<b>Bed-o' Cobs® / Pure-o' Cel®</b>	1/4" Cob & Paper Squares	4PB	115	1.25 cu ft/35 L bag	80
	<b>Enrich-o'Cobs®</b>	1/4" Cob & Paper Rolls	4EB	115	1.25 cu ft/35 L bag	80
	<b>Enrich-n' Pure®</b>	Paper Rolls & Squares	EPB	115	1.25 cu ft/35 L bag	80
NESTING + ENRICHMENT	PRODUCT	DESCRIPTION	PRODUCT CODE	WEIGHT	PACKAGING	PALLET STACK
	<b>Bed-r'Nest®</b> Portion Controlled Nesting	4 gram (kraft)	BRN4SR	15 lb. case	1600/cs (40-ct sleeves)	40
			BRN4BP	30 lb. case	3200/cs (without sleeves)	16
		4 gram (white)	BRN4WSR	15 lb. case	1600/cs (40-ct sleeves)	40
			BRN4WBP	30 lb. case	3200/cs (without sleeves)	16
		8 gram (kraft)	BRN8SR	15 lb. case	800/cs (20-ct sleeves)	40
			BRN8BP	40 lb. case	2200/cs (without sleeves)	16
		8 gram (white)	BRN8WSR	15 lb. case	800/cs (20-ct sleeves)	40
			BRN8WBP	40 lb. case	2200/cs (without sleeves)	16
	<b>Twist-n'Rich™</b>	8" Paper Twists	TNR8	10 lb. case	2500/cs (25 packs of 100)	84
	4" Paper Twists	TNR4	5.7 lb. case	2500/cs (25 packs of 100)	84	
<b>Enrich-n'Nest®</b>	Paper Rolls	EB	-	1.25 cu ft/35 L bag	80	
<b>Crink-I'Nest™</b>	Crinkled Paper (kraft or white)	CNK CNW	22 lb case	Loose in case	15	
<b>Shred-n'Rich™</b>	2 x 4 Enrichment Rectangles	SNR	11.23 lb. case	3000/cs.	32	
SPECIALTY	<b>LabFresh®</b>	Extra Absorbent Litter	LF20	20 lb bag	20 lb bag	77

For more information visit [AndersonsLabBedding.com](http://AndersonsLabBedding.com)

©2020 The Andersons, Inc. All rights reserved. The Andersons logo, Bed-o'Cobs, Bed-r'Nest, Enrich-n'Nest, Enrich-n'Pure, Enrich-o'Cobs, LabFresh, Pel-o'Cobs, and Pure-o'Cel are registered trademarks of The Andersons, Inc. Crink-I'Nest, Twist-n'Rich, Shred-n'Rich, and XTRAFine are trademarks of The Andersons, Inc. 090220

